

Upgrading Your HealthPro® is Easy!

Thank you for upgrading your HealthPro with a V5-Cell filter. This will enable you to upgrade your IQAir HealthPro to a HealthPro Plus in a matter of minutes. As a result, your new HealthPro Plus will filter gaseous pollutants and odors, in addition to allergens, dust and other airborne particles.

1

TM

- A) Ensure that the IQAir device is disconnected from the power supply. Open the two locking arms, ensuring that both are fully extended outwards. (For details, see chapter 6.3 in your Manual.)*
- B) Remove the diffuser and frame 3 containing the IQAir HEPA Filter.
- C) Insert the V5-Cell™ Gas Filter into the empty frame 2 (two). Ensure that the arrows on the filter label point upward. "Top" is indicated on the top of the V5-Cell filter. Press the filter firmly into the frame. (For details, see chapter 6.5 in your Manual.)*
- D) Place frame 3 and the diffuser back on top of the air cleaner.
- E) Push both locking arms inwards simultaneously to close your IQAir system. (For details, see chapter 6.3 in your Manual.)*
- F) Reconnect the power supply.

*User Manual may be downloaded at www.iqair.com/support.


2

Updating the current model setting from HealthPro to HealthPro Plus activates the Filter Life Monitor for the V5-Cell™ filter (Filter 2) on the control panel. This will ensure that the filter life of this new filter will also be properly monitored.


- A) Press the Menu key ten times to reach the Current Model Menu.
- B) Press and hold down the **Enter** key until the cursor appears.
- C) Press the **Left** key to change the model from HealthPro to HealthPro Plus.
- D) Press the **Enter** key to save the HealthPro Plus selection.
- E) Now the control panel will track the filter life of the V5-Cell™ filter. This will be apparent by the green LED lighting up on the air cleaner illustration on the control panel.




Current Model
HealthPro



Current Model
HealthPro



Current Model
HealthPro Plus



Current Model
HealthPro Plus

If you have any questions, please call our Technical Support Team at 888-560-1020 or visit online at www.iqair.com/support.

Thank you for choosing the IQAir® HealthPro® Series.